



# WAVE® Preventative Maintenance Visit

Transgenomic realizes the importance of keeping your laboratory running smoothly and efficiently. Transgenomic WAVE Systems require proper care and maintenance for optimal performance. For laboratories with lower throughput or less defined workload, we offer scheduled preventative maintenance visits.

## SER-PM4500

Preventative Maintenance Visit – WAVE 4500 System

## SER-PM3500

Preventative Maintenance Visit – WAVE 3500 System

## SER-PMMD4000

Preventative Maintenance Visit – WAVE MD 4000 System



These preventative maintenance visits include one day of labor, one preventative maintenance kit and one syringe.

- Replace seals on pump
- Replace flow path filters
- Replace seals on injection valve
- Replace seals on autosampler syringe valve
- Replace syringe
- Check/calibrate oven
- Check lamp energy and wavelength accuracy
- Full clean-up of flow path transfer line and tubing
- Full clean-up of UV detector
- Check leak and pressure sensors
- QC test system to installation standards
- Software updates based on PC configuration and version in use

## Terms and Conditions

- 1) If additional parts are required, they will be ordered separately.
- 2) Travel costs are not included.
- 3) Unless otherwise noted, *Transgenomic Sales Terms and Conditions* apply.

***Discover the Power of Comprehensive Support Services***

**[www.transgenomic.com](http://www.transgenomic.com)**

**Transgenomic Technical Support Center**

**North America**

Phone: (888) 233-9283  
Fax: (402) 452-5453  
Email: [support@transgenomic.com](mailto:support@transgenomic.com)

**Europe**

Phone: +44 141 892 8800  
Fax: +44 141 883 5967  
Email: [eurosupport@transgenomic.co.uk](mailto:eurosupport@transgenomic.co.uk)



**TRANSGENOMIC®**  
*the power of discovery*